

An Interoperability Test Bed for Distributed Healthcare Applications

I-ESA 2014 7th International Conference

Interoperability for Enterprises System and Applications

Albi, France March 24th-28th, 2014

Robert Snelick, National Institute of Standards and Technology (NIST)

March 28th, 2014

Contact: robert.snelick@nist.gov

Overview


- Present an overview of the Interoperability Test Bed
 - Focus on Healthcare Standards
 - Depict a standards based approach to testing
 - Focus on the need for accurate conformance and interoperability requirements
 - Approach is Broadly Applicable
- Describe the Testing Environments
 - Illustrate the use of multiple testing environments
- Identify key aspects of the NIST Testing Infrastructure
 - Depict testing-related capabilities and services
- Present a case study for Patient Identification Cross-referencing
 - Examine the nature of disparate healthcare information technology systems
 - No national ID in the U.S.
 - Demonstrate the operational functionality of the interoperability test bed (ITB) through centrally managed test cases, test flows, test analysis, and test validation reports

NIST – National Institute of Standards and Technology

- U.S. Government Agency under the Department of Commerce
- Information Technology Laboratory (ITL)
 - To promote U.S. innovation and industrial competitiveness by advancing:
 - Measurement science; Standards; and Technology.
 - *through* the research and development in:
 - Information technology; Mathematics; and Statistics.
- Focus on Conformance and Interoperability
 - Specification and Testing of Standards
 - Healthcare Domain
 - Recent mandates within the U.S. to push for interoperability
 - Meaningful Use of EHR Technology (CMS – Center of Medicaid and Medicare)
 - EHR Certification
 - \$40 Billion Dollar Program to promote *meaningful use* of EHR technology
 - Multiple stages (early stages focus is on system capabilities)
 - Later stages to focus on Interoperability—testing platform needed
 - NIST has built numerous conformance testing tools
 - Need for interoperability testing platform

Standards, Conformance, and Testing

- Important of Standards and Testing for Achieving Interoperability
 - “Not just standards but good standards”
 - Clear specification of requirements uses standard conformance constructs
- To achieve interoperability, testing must be an integral part of the process and must be continuous



IHE – Integrating the Healthcare Enterprises

- International organization promoting the use of healthcare standards **integration**
- Integration Profiles (Implementation Guides) that define the workflows in the healthcare domain
 - Data Communication Standards (HL7 V2, V3, CDA, DICOM, etc.)
 - Functional Requirements
- Many Domains
 - Laboratory Orders and Results
 - Infrastructure (e.g., Patient Identification, Document Sharing, etc.)
 - Public Health (Immunization, Syndromic, Reportable Labs, etc.)
 - Etc.
- Testing
 - Pre-connectathon
 - Connectathon
 - Gazelle Test Management Tool

IHE Testing


- Pre-Connectathon
 - Isolated Testing of Vendor's Capabilities
 - Conformance Testing
- Connectathon
 - Peer-to-Peer Testing
 - Interoperability Testing
 - Single Physical Location (U.S., Europe, Japan, etc.)
 - Large Room on closed network
 - Vendor selects domains (Integration Profiles) they wish to test
 - Monitors are used to verify tests (Some automation of tests)
 - Occurs once a year and costly
- IHE Gazelle Test Management Tool
 - Focuses on mostly organizing the test event (not actual testing)
 - Knowledge of systems, integration profiles, creating tests, running tests, evaluating test, evaluating vendors coverage and proficiency of profiles, etc.
 - Some utility to perform some automated testing

Interoperability Test Bed


- Integrated into Gazelle Test Management Tool
- Support Pre-connectathon, Connectathon, and Virtual Connectathon Conformance
- Virtual Connectathon
 - High availability (Internet Based)
 - Leverages Gazelle's utilities to manage test cases
 - Leverages Gazelle's system registration
 - Users can find partners to test with (vendors that implement the same profiles)
 - Publish availability
 - Schedule testing time
 - Integrates testing capabilities and leverages Pre-connectathon conformance test tools and test cases
 - Leverages NIST Testing Infrastructure (used to build standalone conformance test tools)
 - Simulators can be substituted for vendor systems
 - Simulators provide functionality of real system (scaled down reference implementation)
 - Considering stored events of vendors to simulate interoperability testing using simulators
 - Dedicated proxy simplifies test execution

Testing Environments

Isolated System Testing (Conformance Testing)




Peer-to-Peer System Testing (Conformance & Interoperability Testing)




Testing Environments (Replace SUT with Simulator)

Isolated System Testing (Conformance Testing)




Peer-to-Peer System Testing (Conformance & Interoperability Testing)



Interoperability Test Bed

NIST Testing Infrastructure

Services



Test
Artifacts



Validation
Engine




Gazelle Test Management


Test Cases



Integration
Profiles




Configuration




Reports



Test
Execution
Engine




Proxy



Network



Vendor Site A




Vendor **A** Product Under Test

Virtual
Connection




Vendor Site B



Vendor **B** Product Under Test

NIST Testing Infrastructure Core Services




- Provide a portfolio of testing-related capabilities
 - Supported functions include
 - Validation of testing objects
 - Logging of test results
 - Deployment of test agents (Simulators)
 - Generation of test reports
- Preform unique functions
- Well-defined responsibilities and authority
- Multiple Services collaboratively work to support test case execution

NIST Testing Infrastructure


- Centralize health information technology testing resources
- Services-Provide robust support conformance and interoperability testing capability
- Supports a modular test system design
- Flexibility
 - Reusable components
 - Well-defined interfaces
 - Ad-hoc utilization
 - Customizable user defined configurations
- Extensibility
 - Independently ran and operated system components
 - Quick addition of new services
- Scalability
 - Distribution of resource requirements across multiple discrete systems

Isolated System Testing Environment



- Main objective: Conformance testing
 - Conducted with a system under test (SUT) and a testing tool
- Functional behavior assessment
 - Test scenario creates a sequence of transactions to validate requirements
- Test tool
 - Application functionality in normal SUT operational environment
- Upon successful conformance testing, interoperability testing proceeds

Peer-to-Peer System Testing Environment



- Main Objective: Interoperability Testing
 - Conducted with one or more SUTs and a testing tool
- Data Exchange Among a Group of Systems
 - Proxy employed to intercept, log, and rout messages
- Can Leverage Isolated (Conformance) Test Cases
 - Differing test steps, configuration requirements, and assertion assessment
- Interoperability Capability Declarations

CASE STUDY: PATIENT IDENTIFICATION CROSS-REFERENCING


Case Study Overview

- Example case study and workflow that demonstrates typical transactions among disparate healthcare IT systems
- Healthcare systems manage patient data among multiple administrative domains
 - U.S. Patients have many patient identifiers (no national ID)
 - Patients require services in different domains
 - Records may exist across domains
 - Single (master) patient identifier used for document repository retrieval
- Patient identification cross-references is Integration Profile in the IHE infrastructure technical framework
- Case study focuses on using a local patient ID to obtain a master index patient ID to retrieve patient documents from a repository that references the patient information only with the master patient ID
- Resolving patient IDs is used often in practice. One example is querying for a Patient's Immunization Records where the records were submitted to the IIS by multiple sources—getting the correct Patient ID to access the records is the first step of this workflow.
- The Interoperability Ted Bed is design to support workflows for many of the IHE Integration Profiles (e.g., Lab Orders and Results, etc.)


Case Study – Key Components

- Patient identifier and document management systems
 - Patient Identifier Cross-referencing (PIX) Domain
 - Cross-Enterprise Document Sharing (XDS) System
 - Support registering and retrieving documents across enterprise with an administrative domain
- Patient identifier cross-referencing manager (PIX Manager)
 - HL7 V2 data exchange standard
 - Messages enable actor communication
- Actors
 - PIX Source
 - Used for adding and modifying patient demographic data
 - PIX Manager
 - Used for managing and cross-referencing patient identifiers from different domains
 - PIX Consumer
 - Used for querying a PIX Manager for patient identifiers


Patient Identification Cross-Referencing




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued



“Donna Bee” Scenario, continued




“Donna Bee” Scenario, continued




Testing Environments

Isolated System Testing (Conformance Testing)




Peer-to-Peer System Testing (Conformance & Interoperability Testing)




Testing Environments (Replace SUT with Simulator)

Isolated System Testing (Conformance Testing)




Peer-to-Peer System Testing (Conformance & Interoperability Testing)



INTEROPERABILITY TEST BED


NIST Interoperability Test Bed



Logical Operational Divisions / Functions

- **Scheduling:**
 - Coordinates the vendors and aligns interest in test selection, actors, and time availability
- **Configuration:**
 - Records required communication, connection, and addressing for participants
- **Test Setup:**
 - Informs participants of the test case instructions
- **Test Manager:**
 - Controls the overall ITB operation
- **Test Orchestration:**
 - Organizes and coordinates vendor and the test system activities
- **Test Execution:**
 - Workflow management tool; directs test case execution through testing services
- **Status Broadcasting:**
 - Status update broadcast, through served pages, to participants for each performed task
- **Message Handling:**
 - The capturing, storing, analyzing, and forwarding of messages as the progress through the system
- **Singular Analysis:**
 - Each system interaction (message) is assembled and evaluated; a report is rendered and delivered to participants
- **Collective Analysis:**
 - The test analyzer uses the test script for navigating, linking, and analyzing individual validation reports
- **Storing Test Results:**
 - Test results are stored for analysis and auditing purposes


Interoperability Test Bed: Test Flow



Interoperability Test Bed

NIST Testing Infrastructure


Services



Test Artifacts



Validation Engine



Gazelle Test Management


Reports




Test Cases




Configuration



Integration Profiles



Test Execution Engine



Proxy

Network

Vendor Site A

ITB Dashboard

http://

Test Data Specifications



1 Enter data into GUI

Vendor A Product Under Test

Vendor Site B

ITB Dashboard


http://

Vendor A Product Under Test


Interoperability Test Bed

NIST Testing Infrastructure


Services



Test Artifacts



Validation Engine




Gazelle Test Management


Reports




Test Cases




Configuration



Integration Profiles



Test Execution Engine



Proxy

Network

Vendor Site A

Message sent over network

2↑

ITB Dashboard

http://

Test Data Specifications



1 Enter data into GUI

Vendor A Product Under Test

Vendor Site B

ITB Dashboard


http://

Vendor A Product Under Test


Interoperability Test Bed

NIST Testing Infrastructure


Services



Test Artifacts



Validation Engine



Gazelle Test Management

Reports





Test Cases

Configuration

Integration Profiles

Test Execution Engine



Message captured by the Proxy ③

Proxy

Network

Vendor Site A

Message sent over network ②

ITB Dashboard

http://

Test Data Specifications



① Enter data into GUI

Vendor A Product Under Test


Vendor Site B

ITB Dashboard


http://

Vendor A Product Under Test


Interoperability Test Bed




Interoperability Test Bed




Interoperability Test Bed




Interoperability Test Bed




Interoperability Test Bed




Interoperability Test Bed



Interoperability Test Bed



Interoperability Test Bed



INTEROPERABILITY TEST BED DEVELOPMENT

Integrated ITB Functionality into Gazelle

Gazelle - Interoperability Test Bed 2013

Gazelle Registration TF Test Definition Configurations Interoperability Test Bed Administration

Register testing session

Edit Testing Session Information

Year: 2013 Example of value : 2010

Zone: NA Example of values : Europe, NA, Australia

Type: **Interoperability** Example of values : Connection, Showcase

Internet testing: ☐ Testing session is using systems from a not internet testing session

Enable Pre-Connection Testing: ☐ Testing session is using Pre-Connection Testing

Allow Participant Registration: ☒ Checked if Participant Registration shall be enables for that testing session

Test types

Testing session admin

connection HITSP interoperability-test

Copy all Copy Remove Remove All

Copy all Copy Remove Remove All

← Addition of Dropdown to select system defined testing session type.

Test Case Definition Setup

Gazelle - Interoperability Test Bed 2013

Gazelle Registration TF Test Definition Configurations Interoperability Test Bed Administration

Edit Test

Update Test (ITB_PDX_FEED : ITB_PDX_FEED)

back to test list

Test Summary Test Description Test Roles Test Steps

Step Index * 10 Secured ☐ Automatically activate Test Step ☒ Automatically mark Test Step as completed ☒

Initiator Role * PAT_IDENTITY_SRC-PIX

Transaction * ITI-8

Initiator Message Type * ADT_A04

Base Profile for Validation * ADT_A04 - 2.3.1: 1.3.6.1.4.12559.11.1.1.25

Link To Profile ADT_A04 : 1.3.6.1.4.12559.11.1.1.25

Validate Message Using * Please Select...

Validation Context for Initiator Message 2_ORU_R01_2.5.1_ValidationContext.xml Change

Example Message for Initiator MSH|^~^&HL7KITRZ SOFTWARE SERVICES|TEST APP|TEST FACILITY|20090322174532||ADT^A01^ADT_A01|

Option Required

Description Admit/Register a new patient FARNSWORTH^STEVE using these demographics: Name=FARNSWORTH^STEVE, 19781208, Male, Addr=820 JORIE BLVD, CHICAGO, IL 60532, the Patient Identity Source may send ADT^A01 or A04 or A05. This patient is admitted with the Patient ID Assigning Authority value for this Patient Identity Source. (This patient already exists in the Red, Green, and Blue domain)

Responder Role * PIX_MGR-PIX

Responder Message Type * ACK_A04

Base Profile for Validation * Please Select...

Link To Profile ACK_A01 : 1.3.6.1.4.12559.11.1.1.147

Validate Message Using * Gazelle Validation Service

Validation Context for Responder Message VXU_V04_2.5.1_ValidationContext.xml Change

Example Message for Responder MSH|^~^&HL7KITRZ SOFTWARE SERVICES|TEST APP|TEST FACILITY|20090322174532||ACK^A01^ACK_K1|D|2.4|vMSA|AA|v

Associating Testing Artifacts to Test Step

Gazelle - Interoperability Test Bed 2013
tuser1 0
Logout

Gazelle
Registration
TF
Test Definition
Configurations
Interoperability Test Bed
Administration

Test Definition


Test Summary: ITB_PIX_FEED

[back to test list](#)

Test Summary | Test Description | Test Roles | Test Steps





Step Index ▲	Initiator Role ▲	Responder Role ▲	Transaction ▲	Initiator Message Type ▲	Responder Message Type ▲	Option ▲	Description									
10	PAT_IDENTITY_SRC-PIX	PIX_MGR-PIX	ITI-8	ADT_A04	ACK_A04	Required	<p>Admit/Register a new patient FARNSWORTH^STEVE using these demographics: Name=FARNSWORTH^STEVE, 19781208, Male, Addr=820 JORIE BLVD, CHICAGO, IL 60532. the Patient Identity Source may send ADT^A01 or A04 or A05. This patient is admitted with the Patient ID Assigning Authority value for this Patient Identity Source. (This patient already exists in the Red, Green, and Blue domain)</p> <p>Input and Output Contextual Information (0 - 0)</p> <p>Hide</p> <table border="1"> <thead> <tr> <th>Configuration File Type</th> <th>Initiator Scenario Configuration Files</th> <th>Responder Scenario Configuration Files</th> </tr> </thead> <tbody> <tr> <td>Validation Context</td> <td> Download 2_ORU_R01_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre> </td> <td> Download VXU_V04_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre> </td> </tr> <tr> <td>Example Message</td> <td> Download ADT_A01_2.5.1.xml <pre>MSH ^~\& HL7KIT RZ SOFTWARE SERVICES TEST APP TEST FACILITY 20090322174532 ADT^A01^ADT_A01 12345 P 2.5.1 EVN A01 20090322174532 </pre> </td> <td> Download ACK_A01_2.4.xml <pre>MSH ^~\& foo foo foo 200108151718 ACK^A01^ACK 1 D 2.4 v MSA AA^v</pre> </td> </tr> </tbody> </table>	Configuration File Type	Initiator Scenario Configuration Files	Responder Scenario Configuration Files	Validation Context	Download 2_ORU_R01_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre>	Download VXU_V04_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre>	Example Message	Download ADT_A01_2.5.1.xml <pre>MSH ^~\& HL7KIT RZ SOFTWARE SERVICES TEST APP TEST FACILITY 20090322174532 ADT^A01^ADT_A01 12345 P 2.5.1 EVN A01 20090322174532 </pre>	Download ACK_A01_2.4.xml <pre>MSH ^~\& foo foo foo 200108151718 ACK^A01^ACK 1 D 2.4 v MSA AA^v</pre>
Configuration File Type	Initiator Scenario Configuration Files	Responder Scenario Configuration Files														
Validation Context	Download 2_ORU_R01_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre>	Download VXU_V04_2.5.1_ValidationContext.xml <pre><?xml version="1.0" encoding="UTF-8"?> <con:HL7V2MessageValidationContextDefinition xmlns:con="http://www.nist.gov/healthcare/validation/message/hl7v2/co</pre>														
Example Message	Download ADT_A01_2.5.1.xml <pre>MSH ^~\& HL7KIT RZ SOFTWARE SERVICES TEST APP TEST FACILITY 20090322174532 ADT^A01^ADT_A01 12345 P 2.5.1 EVN A01 20090322174532 </pre>	Download ACK_A01_2.4.xml <pre>MSH ^~\& foo foo foo 200108151718 ACK^A01^ACK 1 D 2.4 v MSA AA^v</pre>														
20	PAT_IDENTITY_SRC-PIX	PIX_MGR-PIX	ITI-8	ADT_A01		Required	<p>Admit/Register a new patient SINGLETON^MARION. The Patient Identity Source may send ADT^A01 or A04 or A05. Create patient ID and demographics for Marion Singleton.</p> <p>Input and Output Contextual Information (0 - 0)</p> <p>Configuration Files associated with Test Step</p>									

User Feedback on the Gazelle Dashboard

Test Instance Steps							close	
Step	Trans:	Trans. - Message Type	Opt.	Sending Actor	Receiving Actor	Execution Status	Action	
10		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Waiting		
	Desc:	Admit/Register a new patient FARNSWORTH*STEVE using these demographics: Name=FARNSWORTH*STEVE, 19781208, Male, Addr=820 JORIE BLVD, CHICAGO, IL 60532, the Patient Identity Source may send ADT*A01 or A04 or A05. This patient is admitted with the Patient ID Assigning Authority value for this Patient Identity Source. (This patient already exists in the Red, Green, and Blue domain)					Skip Test Step	
	Logs:	No comment, file or URL  <div>Upload a file (click or drop)</div>					Abort Test Step Pause Test Step	
20		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Inactive		
	Desc:	Admit/Register a new patient SINGLETON*MARION. The Patient Identity Source may send ADT*A01 or A04 or A05. Create patient ID and demographics for Marion Singleton.					Skip Test Step	
	Logs:	No comment, file or URL						

Test Instance Steps

close

Step	Trans:	Trans. - Message Type	Opt.	Sending Actor	Receiving Actor	Execution Status	Action
10		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Skipped	Reactivate Test Step
	Desc:	Admit/Register a new patient FARNSWORTH*STEVE using these demographics: Name=FARNSWORTH*STEVE, 19781208, Male, Addr=820 JORIE BLVD, CHICAGO, IL 60532, the Patient Identity Source may send ADT*A01 or A04 or A05. This patient is admitted with the Patient ID Assigning Authority value for this Patient Identity Source. (This patient already exists in the Red, Green, and Blue domain)					
	Logs:	No comment, file or URL  <div>Upload a file (click or drop)</div>					
20		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Inactive	Activate Test Step Skip Test Step
	Desc:	Admit/Register a new patient SINGLETON*MARION. The Patient Identity Source may send ADT*A01 or A04 or A05. Create patient ID and demographics for Marion Singleton.					
	Logs:	No comment, file or URL  <div>Upload a file (click or drop)</div>					
30		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Inactive	Activate Test Step Skip Test Step
	Desc:	Admit/Register a new patient (ADT*A01 or A04 or A05). This patient has a last name matching the Patient Identity Source's company, and a woman's first name of MARY, eg. patient XYZMedicalCo*MARY. Use this patient name and generate your own patient ID and demographics.					
	Logs:	No comment, file or URL  <div>Upload a file (click or drop)</div>					
40		ITI-8 - ADT_A08	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR	Inactive	Activate Test Step Skip Test Step
	Desc:	Update the patient name (Send ADT*A08 message). Last name stays the same. Change the first name from MARY to MARION. Patient MARION Company_Name exists in PIX Manager, not MARY Company_Name. At this point in the test the PIX Manager contains patients MARION Company_name and MARION SINGLETON.					
	Logs:	No comment, file or URL  <div>Upload a file (click or drop)</div>					
		ITI-8 - ADT_A40	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR		Activate Test Step
	Desc:	SINGLETON*MARION is the maiden name of the patient just registered and updated. Merge SINGLETON*MARION into the last patient updated, eg "XYZMedical*MARION". SINGLETON is subsumed; the patient Marion with the last name that matches your patient name survives (i.e, patient MARION System_name exists in PIX Manager; patient SINGLETON no longer exists).					

Browsing the Proxy for Message Capture

Test Instance Steps

close

Step	Trans:	Trans. - Message Type	Opt.	Sending Actor	Receiving Actor	Execution Status	Action																							
10		ITI-8 - ADT_A01	R	GATEWAY_OR1 (1.0) - PAT_IDENTITY_SRC	GATEWAY_OR2 (1.0) - PAT_IDENTITY_X_REF_MGR																									
	Desc:	Admit/Register a new patient FARNSWORTH^STEVE using these demographics: Name=FARNSWORTH^STEVE, 19781208, Male, Addr=820 JORIE BLVD. CHICAGO, IL 60532. the Patient Identity Source may send ADT^A01 or A04 or A05. This patient is admitted with the Patient ID Assigning Authority value for this Patient Identity Source. (This patient already exists in the Red, Green, and Blue domain)																												
	Logs:	No comment, file or URL																												
	Messages	<div> <div>Upload</div> <div>Message Details</div> <div> <div>Message Details</div> <div> <div>Message ID : 51</div> <div>Date Received : 2/19/14 5:35:03 PM</div> </div> <table> <thead> <tr> <th colspan="2">Initiator</th> <th colspan="2">Proxy</th> <th colspan="2">Responder</th> </tr> </thead> <tbody> <tr> <td>Initiator IP</td> <td>192.168.1.2</td> <td>Proxy IP</td> <td>127.0.0.1</td> <td>Responder IP</td> <td>192.168.1.3</td> </tr> <tr> <td>Initiator Port</td> <td>48558</td> <td>Proxy Port</td> <td>10005</td> <td>Responder Port</td> <td>10056</td> </tr> <tr> <td></td> <td></td> <td>Request Direction</td> <td></td> <td>Message Length</td> <td>465</td> </tr> </tbody> </table> <div>Message Content</div> <div> <pre>MSH ^~\& PAMSimulator IHE PIXManager IHE 20140219043504 ADT^A01^ADT_A01 20140219043504 F 2.3.1 UNICODE UTF-8 EVN 20140219043504 20140219043504 PID DDS- 46929^***DD961.3.6.1.4.1.1.2559.11.1.4.1.26130^PI Turner^Asucena^*****L Foret^*****M 19950612124326 F Race Street^^Gainesville^^Georgia^30501^USA 14^^IHEPAM61.3.6.1.4.1.1.2559.11.1.2.2.5null^AM FV1 16^^IHEPAM61.3.6.1.4.1.1.2559.11.1.2.2.5null^VM V</pre> </div> </div> </div>						Initiator		Proxy		Responder		Initiator IP	192.168.1.2	Proxy IP	127.0.0.1	Responder IP	192.168.1.3	Initiator Port	48558	Proxy Port	10005	Responder Port	10056			Request Direction		Message Length
Initiator		Proxy		Responder																										
Initiator IP	192.168.1.2	Proxy IP	127.0.0.1	Responder IP	192.168.1.3																									
Initiator Port	48558	Proxy Port	10005	Responder Port	10056																									
		Request Direction		Message Length	465																									


Validation Results of a Message Instance

Gazelle - Interoperability Test Bed 2013

tuser10Logout

Gazelle Registration TF Test Definition Configurations Interoperability Test Bed Administration

Validation Result Details



External Validation Report

Disclaimer

The GazelleHL7v2Validator is an experimental system. IHE-Europe assumes no responsibility whatsoever for its use by other parties, and makes no guarantees, expressed or implied, about its quality, reliability, or any other characteristic. We would appreciate acknowledgement if the service is used. Bug tracking service is available at <http://gazelle.ihe.net/jira/browse/HLVAL>

General Information

Validation Date	2014, 02 19 - 10:35 (PM)
Validation Service	Gazelle HL7 Validator (2.2.1)
Message OID	1.3.6.1.4.1.12559.11.1.2.1.5.50438
Profile OID	1.3.6.1.4.12559.11.1.1.24 (Revision: 36704)
Validation Test Status	FAILED

Resources

Find below the list of HL7 tables which have been used to validate the coded values of the message

- 1.3.6.1.4.12559.11.1.3.1.4 (Revision: 36704)

Validation counters

1 error(s)
20 warning(s)
353 report(s)

Validation details

Errors

Type	Required field missing
Location	ADT_A01/PV1[0]/PV1-2(Patient Class)
Description	Element 'Patient Class' is specified as required (R) but not present in the message

Warnings

Type	Length exceeded
Location	ADT_A01/MSH[0]/MSH-3(Sending Application)[0]/HD-1(namespace ID)
Value	DANSimulator

Summary: To achieve interoperability, testing must be an integral part of the standards and development process and testing must be continuous.

